

# Unit Dose Presoak

Enzymatic formula breaks down protein and starch films

## Effective and efficient

- Presoaking flatware results in consistently clean and spotless results with every wash cycle
- Saves time and labor by helping to reduce need for rewash

## Safe to use

- Safer product handling
- Requires no PPE

## Easy on the environment

- Concentrated solid tablets contribute to reduced water, energy, and plastic waste



The tablet's solid chemistry is formulated to save water, energy, and packaging.

**98.8%**

reduction in plastic waste<sup>1</sup>

1. Reduced plastic and cardboard waste enabled by minimizing packaging weight and materials due to solid tablet chemistry. Further detailed information for each product value found is located in disclosure details on Ecolab.com/ReadyDose.

## DIRECTIONS FOR USE:

1. Tear open a packet and drop the tablet into a presoak tub. Wash hands after handling product.
2. Fill up the tub with warm water (110° – 120° F) and allow the tablet to dissolve.
3. Place dirty flatware into the solution and allow to soak.

Unit Dose Presoak	
PACK SIZE	48 ct
DILUTION RATIO	1 tablet/ tub
PICK CODE	6102700
GS CODE	0081900



We thought small so you can think big.  
Ready to learn more about our full tablet cleaning program?  
Visit [cleanwithguestsupply.com](https://cleanwithguestsupply.com)